

INKOLET

LARGE SUPER ABSORBENT INCONTINENCE SHEET



VALUE IN CARING


VENDLOT



PRODUCT INFORMATION

INKOLET is a disposable incontinence sheet with a very large surface and high absorbency.

INKOLET contains fluff pulp, which ensures a good absorbency and rapid spread of liquid. In addition, INKOLET contains superabsorber (SAP), which provides extra absorption and reduces odour. Fluid channels in the sheet ensure that liquid is easily spread into the sheet, which reduces the risk of leaks.

The surface is made of a liquid absorbent non woven material that appears dry and soft against the skin. The surface is anti-bacterial treated and bio-compatible.

The underside is covered with a waterproof PE film. INKOLET is sealed on all four edges. The seal ensures stable and strong edges and ensures that the PE - underside does not come into contact with the patients skin.

INTENDED USE

INKOLET covers a very large surface area with a single sheet. Consequently you achieve sufficient coverage using one large sheet rather than using two or more smaller incontinence sheets. Several incontinence sheets overlapping increases the risk of leaks. Overlaps can also be uncomfortable for the patient and can lead to pressure ulcers.

The large size also makes INKOLET particularly suitable for use with a VENDLET since it is able to extend over the VENDLET-systems bars and thereby follows the patient during turning procedures, protecting sheet and mattress throughout the procedure.

WASHING AND CLEANING



Do not wash

Item no	Product	Size	Absorbency (Rothwell-tested)	Anti-bacterial treatment	Bio-compatible
3601100	INKOLET	90 x 150 cm	3000 ml	Yes	Yes



All rights reserved.
Copyright 2016
Vendlet ApS

Distributor Australia:

HLS Healthcare Pty Ltd
Head Office:
3/56 Norcal Road, Nunawading, VIC 3131
Tel. +61 3 9878 4448

admin@hlshealthcare.com.au
www.hlshealthcare.com.au

Manufacturer:

Vendlet ApS
Egelund 33, 6200 Aabenraa
Denmark
Tel. +45 9631 0050

vendlet@vendlet.dk
www.vendlet.com

