



GH3+ Service module (accessories)

- Time of last safety inspection/service.
- Time of next scheduled safety inspection/service.
- Variable time interval indication of 3-24 months (EU regulation: max. 12-month interval).
- Acoustic notification of next safety inspection/service and visual notification on hand control's LCD display.
- Further information:
 - Indication of time of next replacement of:
 - Lifting strap.
 - Toothed belt.
 - Battery.
- Reset of counter in conjunction with replacement of lifting strap, toothed belt and battery.
- Reset of counter in conjunction with safety inspection/service.
- Disconnection of acoustic signal in the event of no charge.

Changes and updates that can be carried out by technician

- Change of lifting speed (GH3+).
- Upload of new software, regularly released by the supplier and uploaded during annual safety inspection/service.
- Calibration of scale module (GH3+).
- Reset of counter in conjunction with replacement of batteries, lifting strap and toothed belt.
- Disconnection of acoustic signal in the event of no charge.

Arguments

- Ensure that the safety inspection/service is carried out.
- Can monitor when vital parts for operational safety must be replaced, such as lifting strap, battery and toothed belt.
- Increased operational time as new software can be updated on site instead of replacing the PCB as before.
- Adjust the speed to the user.
- Improve maintenance planning.